



## IVTScrip™ mRNA-Human CD244, (Cap 0, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK0106MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK0106MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the CD244 protein. This product can be used in -related researches.

#### Specifications

<b>Modified bases</b>	Pseudo-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_0012917374
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

#### SPECIFICATIONS

<b>Modified bases</b>	Pseudo-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_0012917374
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

## GENE INFORMATION

<b>Alternative Names</b>	h2B4; 2B4
<b>Description</b>	This gene is a protein coding gene. Exhibits MHC class I protein binding activity. Involved in natural killer cell activation involved in immune response; positive regulation of inositol phosphate biosynthetic process; and regulation of cytokine production. Localizes to external side of plasma membrane. Implicated in rheumatoid arthritis.