



IVTScrip™ mRNA-Human CDH11, (Cap 0, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK0235MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK0235MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the CDH11 protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001308392
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001308392
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	OB; OSF-4; ESWS
Description	This gene is a protein coding gene. Predicted to have cadherin binding activity and calcium ion binding activity. Predicted to be involved in several processes, including calcium-dependent cell-cell adhesion via plasma membrane cell adhesion molecules; cell morphogenesis; and cell-cell junction organization. Localizes to extracellular exosome. Implicated in Elsahy-Waters syndrome.