




IVTScrip™ mRNA-Human CD200, (Cap 0, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK0609MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK0609MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the CD200 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_138806.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_138806.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	MRC; MOX1; OX-2; MOX2; CD200
Description	This gene is a protein coding gene. Predicted to have protein binding activity involved in heterotypic cell-cell adhesion. Involved in several processes, including negative regulation of matrix metalloproteinase secretion; positive regulation of arginase activity; and regulation of DNA-binding transcription factor activity. Localizes to axon; cell surface; and neuronal cell body.