




IVTScrip™ mRNA-Human DLX5, (Cap 0, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK0742MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK0742MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the DLX5 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_005221
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_005221
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	SHFM1; SHFM1D
Description	This gene is a protein coding gene. Exhibits transcription regulatory region sequence-specific DNA binding activity. Involved in cell population proliferation and positive regulation of transcription, DN A-templated. Predicted to localize to nuclear chromatin. Implicated in split hand-foot malformation 1 and split hand-foot malformation 1 with sensorineural hearing loss.