



IVTScrip™ mRNA-Human ADARB2, (Cap 0, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK10114MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK10114MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ADARB2 protein. This product can be used in Oogenesis phase fetal germ cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_018702.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_018702.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	RED2; ADAR3
Description	This gene encodes a member of the double-stranded RNA adenosine deaminase family of RNA-editing enzymes and may play a regulatory role in RNA editing. [provided by RefSeq, Jul 2008]