



IVTScrip™ mRNA-Human ABHD13, (Cap 0, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK10262MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK10262MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ABHD13 protein. This product can be used in Leydig cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_032859.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_032859.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	BEM46L1; C13orf6; bA153l24.2
Description	ABHD13 (Abhydrolase Domain Containing 13) is a Protein Coding gene. Diseases associated with ABHD13 include Hypotonia-Cystinuria Syndrome. Gene Ontology (GO) annotations related to this gene include hydrolase activity. An important paralog of this gene is ABHD12B.