



## IVTScrip™ mRNA-Human AC004069.2, (Cap 0, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK10286MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK10286MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the AC004069.2 protein. This product can be used in Natural killer T (NKT) cell-related researches.

#### Specifications

|                |                       |
|----------------|-----------------------|
| Modified bases | 5-Methyl-CTP          |
| 5' Cap         | Cap 0                 |
| Species        | Human                 |
| Applications   | Gene therapy research |
| Format         | Powder                |
| Quantity       | 100 µg                |
| Purification   | Chromatography        |

#### SPECIFICATIONS

|                       |                       |
|-----------------------|-----------------------|
| <b>Modified bases</b> | 5-Methyl-CTP          |
| <b>5' Cap</b>         | Cap 0                 |
| <b>Species</b>        | Human                 |
| <b>Applications</b>   | Gene therapy research |
| <b>Format</b>         | Powder                |
| <b>Quantity</b>       | 100 µg                |
| <b>Purification</b>   | Chromatography        |

## GENE INFORMATION

|                          |                                      |
|--------------------------|--------------------------------------|
| <b>Alternative Names</b> | Nucleoporin 43kDa (NUP43) Pseudogene |
| <b>Description</b>       | ENSG00000251473 is a Pseudogene.     |