



## IVTScrip™ mRNA-Human ACAD10, (Cap 0, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK10304MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK10304MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ACAD10 protein. This product can be used in LGR5+ stem cell-related researches.

#### Specifications

<b>Modified bases</b>	5-Methyl-CTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_001136538.2
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

#### SPECIFICATIONS

<b>Modified bases</b>	5-Methyl-CTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_001136538.2
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

## GENE INFORMATION

### Description

This gene encodes a member of the acyl-CoA dehydrogenase family of enzymes (ACADs), which participate in the beta-oxidation of fatty acids in mitochondria. The encoded enzyme contains a hydrolase domain at the N-terminal portion, a serine/threonine protein kinase catalytic domain in the central region, and a conserved ACAD domain at the C-terminus. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Nov 2008]