



IVTScrip™ mRNA-Human ADAT1, (Cap 0, 5-Methyl-CTP, 30 nt-poly(A))

Cat. No.: GTTS-WK10415MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK10415MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ADAT1 protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001324444.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001324444.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	HADAT1
Description	This gene is a member of the ADAR (adenosine deaminase acting on RNA) family. Using site-specific adenosine modification, proteins encoded by these genes participate in the pre-mRNA editing of nuclear transcripts. The protein encoded by this gene, tRNA-specific adenosine deaminase 1, is responsible for the deamination of adenosine 37 to inosine in eukaryotic tRNA. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jul 2010]