



IVTScrip™ mRNA-Human AC002467.7, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK10885MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK10885MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the AC002467.7 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description LOC101927974 (Uncharacterized LOC101927974) is an RNA Gene, and is affiliated with the lncRNA class.