



IVTScrip™ mRNA-Human ACE2, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK10922MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK10922MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ACE2 protein. This product can be used in Oogenesis phase fetal germ cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001371415.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001371415.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ACEH
Description	The protein encoded by this gene belongs to the angiotensin-converting enzyme family of dipeptidyl carboxydipeptidases and has considerable homology to human angiotensin 1 converting enzyme. This secreted protein catalyzes the cleavage of angiotensin I into angiotensin 1-9, and angiotensin II into the vasodilator angiotensin 1-7. ACE2 is known to be expressed in various human organs, and its organ- and cell-specific expression suggests that it may play a role in the regulation of cardiovascular and renal function, as well as fertility. In addition, the encoded protein is a functional receptor for the spike glycoprotein of the human coronavirus HCoV-NL63 and the human severe acute respiratory syndrome coronaviruses, SARS-CoV and SARS-CoV-2, the latter is the causative agent of coronavirus disease-2019 (COVID-19). Multiple splice variants have been found for this gene and the dACE2 (or MIRb-ACE2) splice variant has been found to be interferon inducible. [provided by RefSeq, Nov 2020]