



## IVTScrip™ mRNA-Human AGBL2, (Cap 1, N1-Methylpseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK11391MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK11391MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the AGBL2 protein. This product can be used in SLC16A7+ cell-related researches.

#### Specifications

|                |                       |
|----------------|-----------------------|
| Modified bases | N1-Methylpseudo-UTP   |
| 5' Cap         | Cap 1                 |
| Species        | Human                 |
| RefSeq         | NM_024783.4           |
| Applications   | Gene therapy research |
| Format         | Powder                |
| Quantity       | 100 µg                |
| Purification   | Chromatography        |

#### SPECIFICATIONS

|                       |                       |
|-----------------------|-----------------------|
| <b>Modified bases</b> | N1-Methylpseudo-UTP   |
| <b>5' Cap</b>         | Cap 1                 |
| <b>Species</b>        | Human                 |
| <b>RefSeq</b>         | NM_024783.4           |
| <b>Applications</b>   | Gene therapy research |
| <b>Format</b>         | Powder                |
| <b>Quantity</b>       | 100 µg                |
| <b>Purification</b>   | Chromatography        |

## GENE INFORMATION

|                          |                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | CCP2                                                                                                                                                                                                                                                                                                        |
| <b>Description</b>       | AGBL2 (AGBL Carboxypeptidase 2) is a Protein Coding gene. Diseases associated with AGBL2 include Cataract 31, Multiple Types and Cardiofaciocutaneous Syndrome 2. Gene Ontology (GO) annotations related to this gene include metallocarboxypeptidase activity. An important paralog of this gene is AGBL3. |