



IVTScrip™ mRNA-Human CD45RO, (Cap 0, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK1156MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1156MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the CD45RO protein. This product can be used in -related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002838
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002838
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180
Description	This gene is a protein coding gene. Exhibits cytoskeletal protein binding activity; protein kinase binding activity; and protein tyrosine phosphatase activity. Involved in several processes, including hematopoietic or lymphoid organ development; positive regulation of cell population proliferation; and regulation of protein kinase activity. Localizes to bleb and external side of plasma membrane. Implicated in hepatitis C; multiple sclerosis; severe combined immunodeficiency; and severe combined immunodeficiency, autosomal recessive, T cell-negative, B cell-positive, Nk cell-positive.