



IVTScrip™ mRNA-Human AGBL4, (Cap 1, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK11692MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK11692MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the AGBL4 protein. This product can be used in Meiotic prophase fetal germ cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001323573.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001323573.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	AGBL4 (AGBL Carboxypeptidase 4) is a Protein Coding gene. Diseases associated with AGBL4 include Cardiofaciocutaneous Syndrome 2 and Cardiofaciocutaneous Syndrome 3. Gene Ontology (GO) annotations related to this gene include tubulin binding and metallocarboxypeptidase activity. An important paralog of this gene is AGBL2.
--------------------	---