



IVTScrip™ mRNA-Human ABHD17A, (Cap 1, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK11766MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK11766MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ABHD17A protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001130111.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001130111.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	C19orf27; FAM108A1
Description	ABHD17A (Abhydrolase Domain Containing 17A, Depalmitoylase) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include hydrolase activity. An important paralog of this gene is ABHD17C.