




IVTScrip™ mRNA-Human ADCYAP1, (Cap 1, 2-Thio-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK11925MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK11925MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ADCYAP1 protein. This product can be used in Astrocyte-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001099733.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001099733.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	PACAP
Description	This gene encodes a secreted proprotein that is further processed into multiple mature peptides. These peptides stimulate adenylate cyclase and increase cyclic adenosine monophosphate (cAMP) levels, resulting in the transcriptional activation of target genes. The products of this gene are key mediators of neuroendocrine stress responses. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2013]