



IVTScrip™ mRNA-Human ADGRG3, (Cap 1, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK12242MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK12242MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ADGRG3 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001308360.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001308360.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ADGRG3 (Adhesion G Protein-Coupled Receptor G3) is a Protein Coding gene. Diseases associated with ADGRG3 include Polymicrogyria, Bilateral Frontoparietal and Epulis. Among its related pathways are Innate Immune System. Gene Ontology (GO) annotations related to this gene include G protein-coupled receptor activity and transmembrane signaling receptor activity. An important paralog of this gene is ADGRG1.
--------------------	---