



IVTScrip™ mRNA-Human ABHD12B, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK12661MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK12661MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ABHD12B protein. This product can be used in Epiblast cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001206673.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001206673.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	BEM46L3; C14orf29; c14_5314
Description	ABHD12B (Abhydrolase Domain Containing 12B) is a Protein Coding gene. Diseases associated with ABHD12B include Chanarin-Dorfman Syndrome. Gene Ontology (GO) annotations related to this gene include hydrolase activity and serine-type peptidase activity. An important paralog of this gene is ABHD12.