




IVTScrip™ mRNA-Human ACTR10, (Cap 1, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK12768MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK12768MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ACTR10 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_018477.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_018477.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	Arp10; Arp11; ACTR11; HARP11
Description	ACTR10 (Actin Related Protein 10) is a Protein Coding gene. Diseases associated with ACTR10 include Brachydactyly, Type A3. Among its related pathways are Innate Immune System and Pathways of neurodegeneration - multiple diseases. An important paralog of this gene is ACTC1.