



IVTScrip™ mRNA-Human ABHD6, (Cap 1, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK12971MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK12971MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the ABHD6 protein. This product can be used in Primitive endoderm cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001320126.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001320126.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ABHD6 (Abhydrolase Domain Containing 6, Acylglycerol Lipase) is a Protein Coding gene. Diseases associated with ABHD6 include Chanarin-Dorfman Syndrome and Retinitis Pigmentosa. Among its related pathways are Signaling by GPCR and Response to elevated platelet cytosolic Ca2+. Gene Ontology (GO) annotations related to this gene include acylglycerol lipase activity. An important paralog of this gene is EPHX3.
--------------------	--