



IVTScrip™ mRNA-Human CD81, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK1317MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1317MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the CD81 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001297649
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001297649
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	TAPA1; TSPAN28
Description	This gene is a protein coding gene. Exhibits cholesterol binding activity; signaling receptor binding activity; and virus receptor activity. Involved in several processes, including positive regulation of cellular metabolic process; positive regulation of immune response; and positive regulation of lymphocyte activation. Localizes to several cellular components, including basal plasma membrane; immunological synapse; and tetraspanin-enriched microdomain. Implicated in common variable immunodeficiency.