



IVTScrip™ mRNA-Human ACBD7, (Cap 1, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK13318MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK13318MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the ACBD7 protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001039844.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001039844.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ACBD7 (Acyl-CoA Binding Domain Containing 7) is a Protein Coding gene. Among its related pathways are Fatty Acyl-CoA Biosynthesis and Regulation of lipid metabolism by Peroxisome proliferator-activated receptor alpha (PPARalpha). Gene Ontology (GO) annotations related to this gene include lipid binding and fatty-acyl-CoA binding. An important paralog of this gene is DBI.
--------------------	---