



## IVTScrip™ mRNA-Human ABCF3, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK13556MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK13556MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ABCF3 protein. This product can be used in SLC16A7+ cell-related researches.

#### Specifications

<b>Modified bases</b>	2-Thio-UTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001351298.2
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

#### SPECIFICATIONS

<b>Modified bases</b>	2-Thio-UTP
<b>5' Cap</b>	Cap 1
<b>Species</b>	Human
<b>RefSeq</b>	NM_001351298.2
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	Chromatography

## GENE INFORMATION

<b>Alternative Names</b>	EST201864
<b>Description</b>	<p>This gene encodes a member of the ATP-binding cassette (ABC) transporter superfamily. ATP-binding cassette proteins transport various molecules across extra- and intracellular membranes. The protein encoded by this gene displays antiviral effect against flaviviruses. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2017]</p>