



IVTScrip™ mRNA-Human TNC, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK1367MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1367MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the TNC protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002160
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_002160
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	GMEM; DFNA56
Description	This gene is a protein coding gene. Exhibits syndecan binding activity. Involved in response to wounding. Localizes to several cellular components, including collagen-containing extracellular matrix; extracellular space; and focal adhesion. Implicated in allergic conjunctivitis; asthma; and autosomal dominant nonsyndromic deafness 56. Biomarker of several diseases, including glioblastoma multiforme; intracranial vasospasm; lung disease (multiple); pleural disease (multiple); and pulmonary hypertension.