




IVTScrip™ mRNA-Human ACY3, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK13681MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK13681MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ACY3 protein. This product can be used in Microglial cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_080658.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_080658.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ACY-3; HCBP1
Description	ACY3 (Aminoacylase 3) is a Protein Coding gene. Diseases associated with ACY3 include Hepatitis C Virus and Canavan Disease. Among its related pathways are Cytochrome P450 - arranged by substrate type and Metabolism. Gene Ontology (GO) annotations related to this gene include identical protein binding and hydrolase activity, acting on carbon-nitrogen (but not peptide) bonds, in linear amides. An important paralog of this gene is ASPA.