




IVTScrip™ mRNA-Human ADGB, (Cap 1, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK13731MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK13731MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ADGB protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_024694.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_024694.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	CAPN16; C6orf103
Description	ADGB (Androglobin) is a Protein Coding gene. Diseases associated with ADGB include Vas Deferens, Congenital Bilateral Aplasia Of. Gene Ontology (GO) annotations related to this gene include iron ion binding and oxygen binding. An important paralog of this gene is CAPN15.