




IVTScrip™ mRNA-Human ANKRD26P3, (Cap 0, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK14612MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK14612MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ANKRD26P3 protein. This product can be used in Meiotic prophase fetal germ cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ANKRD26P3 (Ankyrin Repeat Domain 26 Pseudogene 3) is a Pseudogene.
--------------------	--