



IVTScrip™ mRNA-Human ANKRA2, (Cap 0, 5-Methoxy-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK15195MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK15195MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the ANKRA2 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_023039.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_023039.5
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ANKRA
Description	ANKRA2 (Ankyrin Repeat Family A Member 2) is a Protein Coding gene. Diseases associated with ANKRA2 include Three M Syndrome 1 and Branchiooculofacial Syndrome. Among its related pathways are Signaling events mediated by HDAC Class II. Gene Ontology (GO) annotations related to this gene include histone deacetylase binding and low-density lipoprotein particle binding. An important paralog of this gene is RFXANK.