



## IVTScrip™ mRNA-Human NANOG, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK1533MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK1533MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the NANOG protein. This product can be used in -related researches.

#### Specifications

<b>Modified bases</b>	5-Methyl-CTP & Pseudo-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_001297698.1
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

### SPECIFICATIONS

<b>Modified bases</b>	5-Methyl-CTP & Pseudo-UTP
<b>5' Cap</b>	Cap 0
<b>Species</b>	Human
<b>RefSeq</b>	NM_001297698.1
<b>Applications</b>	Gene therapy research
<b>Format</b>	Powder
<b>Quantity</b>	100 µg
<b>Purification</b>	oligo-dT affinity purification

## GENE INFORMATION

**Description** This gene is a protein coding gene. Exhibits DNA-binding transcription factor activity and RNA polymerase II transcription regulatory region sequence-specific DNA binding activity. Involved in several processes, including endodermal cell fate specification; positive regulation of stem cell proliferation; and somatic stem cell population maintenance. Localizes to nucleolus.