



IVTScrip™ mRNA-Human ANKRD20A9P, (Cap 0, 5-Methyl-CTP, 120 nt-poly(A))

Cat. No.: GTTS-WK15810MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK15810MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP. It encodes the ANKRD20A9P protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases 5-Methyl-CTP

5' Cap Cap 0

Species Human

Applications Gene therapy research

Format Powder

Quantity 100 µg

Purification Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description ANKRD20A9P (Ankyrin Repeat Domain 20 Family Member A9, Pseudogene) is a Pseudogene.