



IVTScrip™ mRNA-Human AK125727, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK16246MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK16246MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the AK125727 protein. This product can be used in AXL+SIGLEC6+ dendritic cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	VASH1-AS1
Description	VASH1-AS1 (VASH1 Antisense RNA 1) is an RNA Gene, and is affiliated with the lncRNA class.