



IVTScrip™ mRNA-Human ANAPC11, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK16368MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK16368MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ANAPC11 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001002244.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001002244.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	APC11; Apc11p; HSPC214
Description	ANAPC11 (Anaphase Promoting Complex Subunit 11) is a Protein Coding gene. Diseases associated with ANAPC11 include Milker's Nodule. Among its related pathways are CDK-mediated phosphorylation and removal of Cdc6 and Activation of cAMP-Dependent PKA. Gene Ontology (GO) annotations related to this gene include ubiquitin-protein transferase activity and ubiquitin-ubiquitin ligase activity.