



IVTScrip™ mRNA-Human ANKRD30A, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK16415MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK16415MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ANKRD30A protein. This product can be used in Luminal epithelial cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_052997.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_052997.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	NY-BR-1
Description	This gene encodes a DNA-binding transcription factor that is uniquely expressed in mammary epithelium and the testis. Altered expression levels have been associated with breast cancer progression. [provided by RefSeq, Nov 2016]