




IVTScrip™ mRNA-Human ANP32B, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK16451MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK16451MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the ANP32B protein. This product can be used in Retinoid acid signaling-responsive fetal germ cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_006401.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_006401.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	APRIL; SSP29; PHAPI2
Description	ANP32B (Acidic Nuclear Phosphoprotein 32 Family Member B) is a Protein Coding gene. Diseases associated with ANP32B include Myopathy, Congenital, Compton-North. Gene Ontology (GO) annotations related to this gene include histone binding and RNA polymerase binding. An important paralog of this gene is ANP32A.