




IVTScrip™ mRNA-Human ANKRD20A9P, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK16710MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK16710MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It encodes the ANKRD20A9P protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ANKRD20A9P (Ankyrin Repeat Domain 20 Family Member A9, Pseudogene) is a Pseudogene.
--------------------	---