



IVTScrip™ mRNA-Human AP1AR, (Cap 0, N1-Methylpseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK16772MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK16772MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & N1-Methylpseudo-UTP. It codes the AP1AR protein. This product can be used in Mitotic arrest phase fetal germ cell-related researches.

Specifications

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001128426.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	N1-Methylpseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001128426.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	2C18; GBAR; C4orf16; PRO0971; gamma-BAR
Description	AP1AR (Adaptor Related Protein Complex 1 Associated Regulatory Protein) is a Protein Coding gene. Diseases associated with AP1AR include Gliosarcoma. Gene Ontology (GO) annotations related to this gene include kinesin binding and AP-1 adaptor complex binding.