



IVTScrip™ mRNA-Human AMIGO2, (Cap 0, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK16952MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK16952MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the AMIGO2 protein. This product can be used in Sertoli cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001143668.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001143668.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ALI1; DEGA; AMIGO-2
Description	AMIGO2 (Adhesion Molecule With Ig Like Domain 2) is a Protein Coding gene. Diseases associated with AMIGO2 include Parkinson Disease 8, Autosomal Dominant and Gastric Adenocarcinoma. An important paralog of this gene is AMIGO1.