



IVTScrip™ mRNA-Human ANAPC15, (Cap 0, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK16969MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK16969MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the ANAPC15 protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001278485.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001278485.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	APC15; HSPC020; C11orf51
Description	ANAPC15 (Anaphase Promoting Complex Subunit 15) is a Protein Coding gene. Diseases associated with ANAPC15 include Deafness, Autosomal Recessive 63 and Autosomal Recessive Non-Syndromic Sensorineural Deafness Type Dfnb. Among its related pathways are CDK-mediated phosphorylation and removal of Cdc6 and Mitotic Metaphase and Anaphase.