



IVTScrip™ mRNA-Human ANKRD13A, (Cap 0, 5-Methoxy-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK17000MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK17000MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methoxy-UTP. It encodes the ANKRD13A protein. This product can be used in Natural killer T (NKT) cell-related researches.

Specifications

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_033121.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methoxy-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_033121.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ANKRD13; NY-REN-25
Description	ANKRD13A (Ankyrin Repeat Domain 13A) is a Protein Coding gene. Diseases associated with ANKRD13A include Anterior Segment Dysgenesis 5 and Atrial Septal Defect 5. An important paralog of this gene is ANKRD13D.