




IVTScrip™ mRNA-Human ANKIB1, (Cap 0, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK17290MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK17290MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ANKIB1 protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_019004.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_019004.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ANKIB1 (Ankyrin Repeat And IBR Domain Containing 1) is a Protein Coding gene. Diseases associated with ANKIB1 include Angiokeratoma Corporis Diffusum With Arteriovenous Fistulas and Cavernous Malformation. An important paralog of this gene is RNF19A.
--------------------	--