




IVTScrip™ mRNA-Human ANKRD10, (Cap 0, 2-Thio-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK17297MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK17297MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 2-Thio-UTP. It encodes the ANKRD10 protein. This product can be used in Exhausted CD8⁺ T cell-related researches.

Specifications

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001286721.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	2-Thio-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001286721.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ANKRD10 (Ankyrin Repeat Domain 10) is a Protein Coding gene. Diseases associated with ANKRD10 include Agraphia. An important paralog of this gene is ANKRD6.
--------------------	--