



IVTScrip™ mRNA-Human CD11b, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK1782MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1782MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the CD11b protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000632.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000632.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	MAC-1; CR3A; MO1A
Description	This gene is a protein coding gene. Exhibits heat shock protein binding activity. Involved in ectodermal cell differentiation; positive regulation of neutrophil degranulation; and positive regulation of superoxide anion generation. Localizes to cell surface; extracellular space; and plasma membrane raft. Implicated in lupus nephritis.