



## IVTScrip™ mRNA-Human ALOX12P2, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK17830MR

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

#### Product overview

This product GTTS-WK17830MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the ALOX12P2 protein. This product can be used in Multilymphoid progenitor cell-related researches.

### Specifications

|                |                           |
|----------------|---------------------------|
| Modified bases | 5-Methyl-CTP & Pseudo-UTP |
| 5' Cap         | Cap 0                     |
| Species        | Human                     |
| Applications   | Gene therapy research     |
| Format         | Powder                    |
| Quantity       | 100 µg                    |
| Purification   | Chromatography            |

### SPECIFICATIONS

|                       |                           |
|-----------------------|---------------------------|
| <b>Modified bases</b> | 5-Methyl-CTP & Pseudo-UTP |
| <b>5' Cap</b>         | Cap 0                     |
| <b>Species</b>        | Human                     |
| <b>Applications</b>   | Gene therapy research     |
| <b>Format</b>         | Powder                    |
| <b>Quantity</b>       | 100 µg                    |
| <b>Purification</b>   | Chromatography            |

## GENE INFORMATION

|                          |                                                                       |
|--------------------------|-----------------------------------------------------------------------|
| <b>Alternative Names</b> | ALOX12E                                                               |
| <b>Description</b>       | ALOX12P2 (Arachidonate 12-Lipoxygenase Pseudogene 2) is a Pseudogene. |