



IVTScrip™ mRNA-Human ANKRD28, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK17914MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK17914MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the ANKRD28 protein. This product can be used in Meiotic prophase fetal germ cell-related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001195098.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001195098.1
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	PITK; FAP79; CFAP79; PPP1R65
Description	ANKRD28 (Ankyrin Repeat Domain 28) is a Protein Coding gene. Diseases associated with ANKRD28 include Urethra Clear Cell Adenocarcinoma and Urethra Adenocarcinoma. Among its related pathways are Vesicle-mediated transport and Metabolism of proteins. An important paralog of this gene is ANKRD44.