



IVTScrip™ mRNA-Human ANKRD44, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK17926MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK17926MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the ANKRD44 protein. This product can be used in Lymphoid-primed multipotent progenitor cell-related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001195144.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_001195144.2
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	PP6-ARS-B
Description	ANKRD44 (Ankyrin Repeat Domain 44) is a Protein Coding gene. Diseases associated with ANKRD44 include Blepharoconjunctivitis and Glass Syndrome. An important paralog of this gene is ANKRD28.