




IVTScrip™ mRNA-Human CK14, (Cap 0, 5-Methyl-CTP & Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK1796MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1796MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 0 & 5-Methyl-CTP & Pseudo-UTP. It encodes the CK14 protein. This product can be used in -related researches.

Specifications

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000526
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	5-Methyl-CTP & Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_000526
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Alternative Names	EBS3; EBS4
Description	This gene is a protein coding gene. Exhibits keratin filament binding activity. Involved in aging; hair cycle; and intermediate filament bundle assembly. Localizes to cytoplasm; keratin filament; and nucleus. Implicated in Naegeli-Franceschetti-Jadassohn syndrome; dermatopathia pigmentosa reticularis; and epidermolysis bullosa simplex (multiple).