



IVTScrip™ mRNA-Human AHSA2, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK18027MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK18027MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the AHSA2 protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_152392.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_152392.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	AHA1; Hch1; AHSA2P
Description	AHSA2P (Activator Of HSP90 ATPase Homolog 2, Pseudogene) is a Pseudogene. Diseases associated with AHSA2P include Chromosome 2P16.1-P15 Deletion Syndrome and Palmoplantar Keratoderma And Congenital Alopecia 2.