



IVTScrip™ mRNA-Human ALKBH4, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK18125MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK18125MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ALKBH4 protein. This product can be used in SLC16A7+ cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_017621.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_017621.4
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Alternative Names	ABH4
Description	ALKBH4 (AlkB Homolog 4, Lysine Demethylase) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include actin binding and oxidoreductase activity, acting on paired donors, with incorporation or reduction of molecular oxygen, 2-oxoglutarate as one donor, and incorporation of one atom each of oxygen into both donors.