



IVTScrip™ mRNA-Human ANKRD65, (Cap 1, Pseudo-UTP, 120 nt-poly(A))

Cat. No.: GTTS-WK18234MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK18234MR is a type of mRNA having 120 nt poly(A) tail and modified with Cap 1 & Pseudo-UTP. It encodes the ANKRD65 protein. This product can be used in Ciliated cell-related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001145210.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 1
Species	Human
RefSeq	NM_001145210.3
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	Chromatography

GENE INFORMATION

Description	ANKRD65 (Ankyrin Repeat Domain 65) is a Protein Coding gene. Diseases associated with ANKRD65 include Brachial Plexus Neuritis and Osteosclerotic Myeloma. An important paralog of this gene is TANC1.
--------------------	--