




IVTScrip™ mRNA-Human SUSD2, (Cap 0, Pseudo-UTP, 30 nt-poly(A))

Cat. No.: GTTS-WK1828MR

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product overview

This product GTTS-WK1828MR is a type of mRNA having 30 nt poly(A) tail and modified with Cap 0 & Pseudo-UTP. It encodes the SUSD2 protein. This product can be used in -related researches.

Specifications

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_019601
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

SPECIFICATIONS

Modified bases	Pseudo-UTP
5' Cap	Cap 0
Species	Human
RefSeq	NM_019601
Applications	Gene therapy research
Format	Powder
Quantity	100 µg
Purification	oligo-dT affinity purification

GENE INFORMATION

Description	This gene is a protein coding gene. Involved in negative regulation of cell cycle G1/S phase transition and negative regulation of cell division. Localizes to plasma membrane.
--------------------	---